

DICOM Viewer Quick Start Job Aid

Study Search and Filtering

- Search and filter for specific patients and studies
- Select a tab to use a pre-defined worklist (worklists vary by site)
- Click the door icon to logout
- Click on the column headers to sort the worklist

Note: If a patient is marked as confidential and you do not have user privileges to view confidential patients, then patients marked as confidential are not included in your search results or in Study Lists.

The screenshot shows the 'Study Search' interface. On the left is a 'Filter' sidebar (A) with fields for Patient Name, Patient ID, Accession #, Date Of Birth, Referring Physician, Exam Status, Study Description, Location, Referring Service, and Modality. Below the filter is an 'Exam Date' section with radio buttons for 'Today', 'Yesterday', and 'Last 7 days', and a 'Date Range' section. At the bottom of the filter are 'Starts' and 'Ends' date pickers and an 'Include' checkbox for 'Zero Image Studies'. A 'Reset' button (G) is at the bottom left of the filter. The main area has tabs for 'All Studies-OPACS', 'CT-ALL', 'OPACS-MR', and 'OPACS-CR' (B). Below the tabs is a table with columns: Name, NHS #, Patient ID, Date Of Birth, Description, Domain, Accession #, Referring Physician, Exam Status, Exam Date, and Modality (D). A table row (E) is highlighted. A 'Search' button (F) is at the bottom center. A door icon (C) is in the top right corner.


- Single-click anywhere in the study row to view the study
- Click Search to view all studies that meet the filtering criteria
- Click reset to clear all filters or start new search

For questions or help with Seattle Radiology PACS,
please contact 206.233.7700 or helpdesk@Radiax.com

DICOM Viewer Quick Start Job Aid


ZFP DICOM Viewer

- A. Access to other studies in the patient history
- B. Reports and Exam Notes open automatically for the Primary study
- C. Click **X** to close the Reports/Exam Notes tab and view images in the full viewer



- D. Key Images for the Primary study open automatically
- E. Click and drag series thumbnail to a viewport

- A. Click to return to the Search page
- B. Change viewport layout, access measurements and image manipulation tools, and launch 3D tools
- C. Scroll using scroll wheel or left-mouse click and drag
- D. Click **?** to open a drop-down menu: User Manual option opens the ZFP online help; Gesture Help option opens the Gesture Help page



- E. Right-click viewport to access further tools and options
- F. CT W/L presets in right-click menu
- G. Double-click an image to 1-up
- H. Right-mouse click and drag to change W/L; press and drag middle mouse to zoom an image